



# In vitro evaluation of daily vs weekly denture cleansing regimens

David Bradshaw<sup>1</sup>, Munisha Mangal<sup>1</sup>, Tracy Young<sup>2</sup>, Jason Brown<sup>2</sup>, Gordon Ramage<sup>2</sup>

<sup>1</sup>GSK Consumer Healthcare, Weybridge, UK

<sup>2</sup>Glasgow Dental Hospital & School, Glasgow, Scotland, UK

## Objectives


To compare the efficacy of daily denture cleanser (DDC) treatment with brushing using standard toothpastes or no treatment in removing and killing *in vitro* denture biofilms.

## Methods

- Biofilms were developed for 7 days on on PMMA discs<sup>(1)</sup>
- Comprised *S.mitis*, *S.intermedius*, *S.oralis*, *C.albicans*, *A.naeslundii*, *V.dispar*, *F.nucleatum* and *F.nucleatum ssp. vincentii*, *P.intermedia*, *P.gingivalis* and *A.actinomycetemcomitans*
- Biofilms were treated by
  - daily brushing and soaking with 3-minute exposure to denture cleanser (DDC)
  - brushing with wetted toothbrush, denture cleanser treatment on D7 only (BR);
  - daily brushing & soaking in toothpaste 1 slurry (25%w/w in hard water) – (TP1)
  - daily brushing and soaking in toothpaste 2 slurry (25%w/w in hard water) (TP2)
  - no treatment (NT) – negative control.

Biofilm viability was then measured by viable counting (VC, total aerobes, anaerobes and yeasts), qPCR (colony forming equivalents, CFE) and visualised using confocal laser scanning microscopy (CLSM).

Figure 1: Experimental Treatments – Daily, Weekly, Toothpaste 1, Toothpaste 2, Control



## Results

Treatments were assessed by comparing the geometric means of VC or CFE.

- NT biofilms contained a geometric mean  $1.99 \times 10^9$  CFU/ml.
- DDC treatment reduced this to  $8.24 \times 10^5$  CFU/ml (2414-fold reduction).
- BR treatment reduced VC by 2.3-fold compared to NT.
- TP1 treatment reduced biofilm CFU by 10.1-fold vs NT, whereas TP2 reduced biofilms by 12.1-fold vs NT.
- Comparing DDC treatment to other treatments showed that this produced:
  - a 1034-fold reduction compared to BR
  - a 240-fold reduction compared to TP1
  - a 200-fold reduction vs TP2.
- Geometric means of total CFU of DDC treated biofilms were significantly better compared with all other treatments ( $P < 0.00041$  or less for each comparison, t-test).
- CFE analysis showed similar, significant benefits for DDC compared with other treatments.

## Results

Figure 2: Daily Treatments CFU





Figure 3: Fold Change in CFU vs Untreated Control



Figure 4: Representative CLSM Images of Treated Biofilms



## Discussion

- Daily denture cleanser treatment gave a significantly greater reduction of biofilms than any other treatment or untreated control.
- Future studies could investigate the effects of toothpaste slurry effects on denture materials which although widely used by consumers can be potentially damaging for dentures<sup>(2)</sup>.

## Conclusions

- Clear numerical differences (by CFU/CFE and visual differences (by CLSM) were observed between treated and untreated biofilms.
- Denture cleanser applied with brushing is a significantly more effective strategy for managing denture biofilms in comparison with no treatment or with intermittent treatment.
- Daily modality of treatment was more effective than weekly.
- Brushing with toothpaste was only moderately effective.

## References

- O'Donnell LE, Robertson D, Nile CJ, Cross LJ, Riggio M, Sherriff A, Bradshaw D, Lambert M, Malcolm J, & MJ, Zaura E, Crielaard W, Brandt BW, Ramage G (2015). The Oral Microbiome of Denture Wearers Is Influenced by Levels of Natural Dentition. *PLoS One* 10(9):e0137717.
  - Axe AS, Varghese R, Bosma M, Kitson N, Bradshaw DJ. (2016). Dental health professional recommendation and consumer habits in denture cleansing. *Journal of Prosthetic Dentistry* 115:183-188.
- This study was funded by GSK Consumer Healthcare. Bradshaw & Mangal are employees of GSK Consumer Healthcare.